

(201) 825-9500 / (800) 631-3826

Vertical Drop CII Sharps Shield



Capintec offer the most complete line of Sharps Shields and containers available for Nuclear Medicine. These Sharps leaded container cabinets offer a convenient way to store spent syringes in a safe manner. Two sizes are shown: One is designed for containment of one medium-size Chimney-Top Monoject® Sharps container and the other holds two small or one large Monoject Sharps Container.

The Vertical Drop CII Sharps Shield comes with a Key-Lock safety feature not found elsewhere. The hinged top is easily opened to remove the Sharps container for replacement, requiring no removal of a heavy lid. Shielded with 1/8" lead, this unit can be specified with 1/4" shielding.

The Vertical Drop CII Sharps is designed to use a Chimney-Top Sharps Container which protects the user's fingers in use.. Convenient needle notches facilitate safe and rapid removal of needles from tube holders.

Features

- All units have Key-Lock tops to meet OSHA requirements
- Table Top or Flush Mount in work counter
- Single or double stacked models
- Lead thickness 1/8" or 1/4"
- Hinged top for easy opening and removal of sharps

Specifications

- Standard Model: Shielding: 1/8" lead with 4 pi shielding; Size: 7.25" W x 11" L x 11" H; Weight: 54lb.
- Tall Model: Shielding: 1/8" lead with 4 pi shielding; Size: 7.25" W x 11" L x 18" H; Weight: 76lb.

Vertical Drop CII Sharps Shield (Standard)-

Item #0660-0018

Vertical Drop CII Sharps Shield (Tall)- Item #0660-0019

Monoject Sharps Medium (For Standard) Box of 20..... 0660-1811

Monoject Sharps Large (For Tall) Box of 10 0660-1812

Monoject Sharps Small (For Tall) Box of 40 0660-1813

Horizontal Drop CII Sharps Shield



These Sharps leaded containers are available for Nuclear Medicine needs. Shown are: one designed for containment of one medium size Horizontal-Entry Monoject® Sharps Container or one medium Gator Sharps container; the other holds one tall Gator Sharps container.

The Horizontal Drop CII Sharps Shield comes with a Key-Lock safety feature, and hinged top and easy opening for replacement. Shielded with 1/8" lead, this unit can be specified with 1/4" shielding.

The Horizontal Drop CII Sharps is designed to use a Horizontal-Entry Sharps Container with greater filling capacity since the syringes lie flat.

Features

- All units have Key-Lock tops to meet OSHA requirements
- Lead thickness 1/8" or 1/4"
- Hinged top for easy opening and removal of sharps

Specifications

- Standard Model: Shielding: 1/8" lead with 4 pi shielding ; Size: 7.25" W x 11" L x 11" H; Weight: 52lb.

Horizontal Drop CII Sharps Shield (Standard) 0660-0038

Davon Gator Sharps Medium (For Standard) case of 30..... 0660-1832